



### Surface disinfectant for use in vaporiser

- Residual antimicrobial effect for up to 5 days
- effective against Coronavirus
- pleasant odour



#### Intended use

**Bacoban<sup>DL</sup> 3%** is an alcohol-free ready-to-use disinfecting solution with residual antimicrobial effect for use in our vaporiser.

#### Only for use in combination with our vaporiser.

Remove the bottle cap and screw the bottle on the inner thread of the vaporiser diffuser attachment carefully and tightly. For further information please follow the instructions in the manual.

#### Efficacy

Testing according to EN14476 (dirty conditions)	time of exposure
Virucidal activity against enveloped viruses (incl. <b>Coronavirus</b> , HIV, HBV, HCV, Influenza, Herpes, Ebolavirus...)	5 min
Rotavirus	5 min
Norovirus	30 min
Human papillomavirus (HPV)	60 min
Adenovirus	240 min

#### EN guidelines

EN 1040 (bactericidal)	5 min
EN 1275 (yeasticidal)	5 min
EN 1276 (bactericidal)	5 min
EN 1650 (yeasticidal)	5 min
EN 13697 (bactericidal)	15 min
EN 13727 (bactericidal)	5 min
EN 13624 (yeasticidal, fungicidal)	15 min
EN 14561 (bactericidal)	5 min
EN 14562 (yeasticidal)	60 min
EN 16615 (bactericidal, yeasticidal) (4 field test)	5 min

#### VAH list (tested according to VAH-guidelines 2015)

bactericidal and yeasticidal	5 min
------------------------------	-------

**Residual antimicrobial effect for up to 5 days** (tested according to ASTM E 2180).

#### Please note:

Do not stay in the same room whilst disinfection is being carried out with the vaporiser. Keep windows and doors closed.

#### Ingredients

100 g of solution containing: Benzalkonium chloride 0,78 g, sodium pyrrhione 0,075 g, polycondensates, perfume substances, purified water.

Does not contain formaldehyde. Does not contain phenols.

#### Storage

Store **Bacoban<sup>DL</sup> 3%** at +5 to +25°C. Protect from frost and direct sunlight.

#### Type of delivery

Product	Packaging unit	Item number
<b>Bacoban<sup>DL</sup> 3%</b> Bottle	6 x 1000 ml	BAC1DL3
<b>Bacoban<sup>DL</sup> 3%</b> Canister	1 x 20 l	BAC20DL3

#### Test reports

Our products are tested extensively. We are pleased to provide you with our microbiological reports on request.

#### Safety notice

Safety data sheet available on request.

**Important!** The residual antimicrobial effect of **Bacoban<sup>DL</sup> 3%** is no excuse for not carrying out necessary disinfection measures in cases of visible and/or questionable contamination!

**Use disinfectants safely. Always read label and product information before use.**

Biocide Reg. No: N-37828 (D)