



## Ready-to-use solution for alcohol-free surface disinfection and cleaning

- Residual antimicrobial effect for up to 10 days
- effective against Coronavirus
- excellent skin tolerance
- high cleaning effect
- practical applicability tested in the 4-field-test



### Intended use

**Bacoban<sup>DL</sup>** is a ready-to-use, waterbased solution for surface disinfection and cleaning with residual, antimicrobial effect.

### Directions for use

Apply **Bacoban<sup>DL</sup>** onto the surface and spread evenly with a cloth respectively take a ready-to-use wipe from the package and spread the liquid evenly on the surface. Then let dry. Don't rinse or rub dry.

### Efficacy

Testing according to EN14476 and guideline DVV / RKI	time of exposure
Virucidal activity against enveloped viruses (incl. <b>Coronavirus</b> , HIV, HBV, HCV, Influenza, Herpes, Ebolavirus...)	5 min
Rotavirus	5 min
Norovirus	30 min
Human papillomavirus (HPV)	60 min
<b>VAH list (tested according to VAH-guidelines 2015)</b>	
bactericidal and yeasticidal	5 min
<b>EN guidelines</b>	
EN 1040 (bactericidal)	5 min
EN 1275 (yeasticidal)	5 min
EN 1276 (bactericidal)	5 min
EN 13727 (bactericidal)	5 min
EN 13624 (yeasticidal, fungicidal)	5 min
EN 1650 (yeasticidal)	5 min
EN 13697 (bactericidal)	15 min
EN 14561 (bactericidal)	60 min
EN 14562 (yeasticidal)	60 min
EN 16615 (bactericidal, yeasticidal) (4 field test)	5 min

**Residual antimicrobial effect for up to 10 days** (tested according to ASTM E 2180).  
Dermatologically tested (epicutaneous test with 30 test persons): excellent skin tolerance.

### Ingredients

100 g of solution containing: Benzalkonium chloride 0,26 g, sodium pyrithione 0,025 g, polycondensates, perfume substances, purified water.

Ingredients according to regulation (EC) No 648/2004: less than 5 % cationic surfactants, disinfectants, perfumes. Does not contain formaldehyde. Does not contain phenols.

### Storage

Store **Bacoban<sup>DL</sup>** at +5 to +25°C. Protect from frost and direct sunlight.

### Type of delivery

Product	Packaging unit	Item number
<b>Bacoban<sup>DL</sup></b> Bottle	12 x 500 ml	BAC500DL
<b>Bacoban<sup>DL</sup></b> Wipes	16 x 50 pcs	BACDLTUE
<b>Bacoban<sup>DL</sup></b> Canister	1 x 5 l	BAC5000DL1

### Test reports

Our products are tested extensively. We are pleased to provide you with our microbiological and dermatological reports on request.

### Safety notice

Safety data sheet available on request.

**Important!** The residual antimicrobial effect of **Bacoban<sup>DL</sup>** is no excuse for not carrying out necessary disinfection measures in cases of visible and/or questionable contamination!

**Use disinfectants safely. Always read label and product information before use.**

Biocide Reg. No: N-34071 (D), 262/17/L (Lu)