



Alcohol-free Disinfection and Cleaning

- Residual antimicrobial effect for up to 10 days
- VAH-listed
- effective against Coronavirus
- excellent skin tolerance
- high cleaning effect
- practical applicability tested in the 4-field-test



VAH-listed



dermatest-guarantee.de



Intended use

Bacoban^{WB} is a waterbased concentrate for combined disinfection and cleaning of inventory and floors with residual antimicrobial effect.

Directions for use

Prepare a ready-to-use solution from **Bacoban^{WB}** according to the mixing table and moisten the surface sufficiently. Don't rinse or rub dry. Lifetime of freshly prepared working solution: 8 hours

Final volume	Mixing table			
	0,5 %	1 %	2 %	3 %
1 Litre	5 ml	10 ml	20 ml	30 ml
2 Litre	1 ml	20 ml	40 ml	60 ml
4 Litre	20 ml	40 ml	80 ml	120 ml
8 Litre	40 ml	80 ml	160 ml	240 ml
10 Litre	50 ml	100 ml	200 ml	300 ml

Residual antimicrobial effect for up to 10 days (tested according to ASTM E 2180).
Dermatologically tested (*at 1 %): excellent skin tolerance.

Storage

Store **Bacoban^{WB}** at +5 to +25°C. Protect from frost and direct sunlight.

Type of delivery

Product	Packaging unit	Item number
Bacoban^{WB} Bottle	6 x 1 l	BACWB1D
Bacoban^{WB} Canister	1 x 5 l	BACWB5

Test reports

Our products are tested extensively. We are pleased to provide you with our microbiological and dermatological reports on request.

Safety notice

Safety data sheet available on request.

Important! The residual antimicrobial effect of **Bacoban^{WB}** is no excuse for not carrying out necessary disinfection measures in cases of visible and/or questionable contamination!

Ingredients

100 g of solution containing: Benzalkonium chloride 26 g, sodium pyrithione 2,5 g, polycondensates, perfume substances, purified water.

Ingredients according to regulation (EC) No 648/2004: 15 % or over but less than 30 % cationic surfactants, disinfectants, perfumes (Linalool, Limone, Citral, Geraniol, Citronellol). Does not contain formaldehyde. Does not contain phenols.

Efficacy

	dilution	time of exposure
Testing according to EN14476 and guideline DVV / RKI		
Virucidal activity against enveloped viruses (incl. Coronavirus , HIV, HBV, HCV, Influenza, Herpes, Ebolavirus...)	1,00%	5 min
Adenovirus	2,00%	240 min
Human papillomavirus (HPV)	1,00%	60 min
Norovirus	1,00%	30 min
Rotavirus	0,10%	5 min
VAH list (tested according to VAH-guidelines 2015)		
bactericidal und yeasticidal	1,00%	5 min
EN guidelines		
EN 1040 (bactericidal)	0,25%	5 min
EN 1275 (yeasticidal)	0,25%	5 min
EN 1276 (bactericidal)	0,75%	5 min
	0,50%	15 min
EN 13727 (bactericidal)	1,00%	5 min
EN 13624 (yeasticidal, fungicidal)	1,00%	5 min
EN 1650 (yeasticidal)	0,50%	5 min
	0,25%	15 min
EN 13697 (bactericidal)	0,50%	15 min
EN 14561 (bactericidal)	3,00%	5 min
	0,50%	60 min
EN 14562 (yeasticidal)	0,50%	60 min
EN 16615 (bactericidal, yeasticidal) (4 field test)	1,00%	5 min

Use disinfectants safely. Always read label and product information before use.

Biocide Reg. No: N-28795 (D). 263/17/L (Lu)